

Quick Guide - Cell Viability Trypan Blue Staining


Contents

- 1x 20 mL Trypan blue 0.02%
- 1x 20 mL Phosphate Buffered Saline w/o Ca^{2+} & Mg^{2+} (DPBS⁻)

Protocol for 96-well Half Area Plates

Example for a 1:20 dilution and 4 wells per sample:

1. Pipet 90 μL DPBS⁻ into a e.g. 1.5 mL centrifuge tube
2. Add 10 μL cell suspension
3. Add 100 μL Trypan blue staining solution and mix well but gently
4. Pipet 40 μL into 4 wells of a 96-well half area plate
5. Repeat steps 1. – 4. for every sample
6. Centrifuge plate 1 min at 30 x g
7. Measure plate



Quick Guide
QG-B004-XVII-01

Expected Cell Density	Final Dilution	PBS ⁻ [μL]	Cell Sample [μL]	Trypan Blue [μL]
$< 0.1 \times 10^6$	1:2	-	100	100
$0.1 - 0.5 \times 10^6$	1:5	60	40	100
$0.5 - 1 \times 10^6$	1:10	80	20	100
$1 - 3 \times 10^6$	1:20	90	10	100
$3 - 6 \times 10^6$	1:40	95	5	100
$6 - 10 \times 10^6$	1:80	97.5	2.5	100

