


Suspension Cell Count (1F)

General Purpose

The Suspension Cell Count (1F) analysis algorithm operates with a brightfield and one fluorescence channel. The image analysis counts all cells in the brightfield image and searches in a second step for a fluorescence signal in the area of the detected cells. It can be used as an endpoint determination e.g. with a viable (e.g. Calcein-AM) or dead (e.g. Propidium Iodide) cell marker to check the culture condition. Furthermore determination of transfection efficiency with marker-co-transfection is another application.



Short Note
SN-B121-XVII-04

Result Table

• TC-BF [#]	Total number of cells detected in the brightfield channel
• TC-1F [#]	Total number of cells in the 1st fluorescence channel
• Living Cells [#]	Quantity of living cells
• Viability	Percentage of viable cells in your sample
• CD-1F [#/mL]	Number of cells labeled with fluorescence 1 per sample volume
• 1F/BF [%]	Fluorescence ratio: $(TC-1F/TC-BF) * 100$
• CD [#/mL]	Cell density = Total number of cells per sample volume
• Sample Volume [μ L]	'Undiluted' volume of your sample
• Avg Cell Size [μ m ²]	Average cell size per well
• Avg Fluo CH1 Intensity BC	Average fluorescence intensity in channel 1 over background
• Sum of Cell Sizes [μ m ²]	Sum of cell sizes

Example

This examples show suspension cells with a green fluorescent CD10 marker.

Marked green:

Detected in brightfield only

Marked orange:

Detected in BF AND fluorescence channel




Fig. 1: a) Overlay of brightfield and fluorescence channel with image analysis. b) Brightfield image of suspension cells. c) Fluorescence channel (FITC-labeled CD10-antibody).